



SYSTEMATIC HYALURONIC ACID GELS



amalian

Dosage recommendation for physicians

MONOPHASIC,
FULLY HYDRATED HYALURONIC
ACID FILLERS
"SF"



amalian SF 16 fine



amalian SF 20 medium



amalian SF 24 advanced



amalian SF 14 natural lips

amalian SF 16 fine



INDICATIONS:

- 16mg/ml, soft cross-linked hyaluronic acid
- fine lines and wrinkles
- glabella
- laugh lines
- forehead lines
- easy nasolabial - and marionette lines
- vermillion border
- augmentation of the hands

Important information:

This gentle hyaluronic acid filler for easy augmentation is ideal for fine wrinkles, vermillion border, forehead lines, crow's feet, light nasolabial and marionette lines.

This fully hydrated gel is evenly distributed in the tissues, can be easily and gently injected and achieves a fresh, young finite appearance while preserving the naturalness.

Like all SF-products is also well tolerated amalian SF 16 fine. Side effects such as swelling are virtually eliminated. The injection should be in the upper to middle dermis.

amalian SF 16 fine is also ideally suited to perform a post-treatment after injection with the stronger product amalian SF 24 advanced. Thus, the volume of building effect is amplified and extended the residence time in the skin.

Needle size: 27G

fine forehead lines


glabella

fine nasolabial lines

laugh lines

fine marionette lines

vermillion border



amalian SF 20 medium



INDICATIONS:

- 20mg/ml, soft cross-linked hyaluronic acid
- nasolabial folds
- marionette lines
- light scars
- increased perioral wrinkles
- neck wrinkles
- décolleté


Important information:

This hyaluronic acid gel is an "allround" filler for medium wrinkles and is injected into the upper to middle layer of skin. Thanks to the excellent flow properties of the gel injections amalian SF 20 medium are very gentle and atraumatic.

The result is a natural lifting effect.

This fully hydrated gel is suitable for almost any application, with the exception of the glabella. The soft texture of the gel homogeneously distributed in the skin, which fills wrinkles on evenly and immediately produces a visible effect. amalian SF 20 medium can also be used with very good success for middle lip augmentation.

Needle size: 27G



light scars

nasolabial folds

marionette lines

neck wrinkles

décolleté

amalian SF 24 advanced



INDICATIONS:

- 24mg/ml, soft cross-linked hyaluronic acid
- medium to very deep wrinkles
- correction of facial imperfections
- development of facial defects
- stronger scars
- nose corrections
- cheek augmentation
- marionette wrinkles
- nasolabial wrinkles

Important information:

amalian SF 24 is a monophasic advanced fully hydrated filler, volumizer, efficient and excellent tolerability. amalian SF 24 is ideally suited for advanced injections deeper folds, to scars correction and modeling of facial contours. Skin irregularities, for example acne scars are smoothed, nasolabial wrinkles are lifted visible. This filler is also **very** suitable for long-lasting hand-augmentation, which is best done in these cases with a blunt cannula. Also for a cheek augmentation it is recommended to use a blunt cannula. amalian SF 24 advanced is injected into the mid to deep skin layer.

Needle size: 27G

correction of
facial imperfections

development of
facial defects


cheek augmentation

marionette wrinkles

stronger scars

nose corrections

nasolabial wrinkles



amalian SF 14 natural lips



INDICATIONS:

- 14mg/ml, soft cross-linked hyaluronic acid
- light lip augmentation

Important information:


amalian SF 14 natural lips is a low concentrated filler for light lip augmentation. It reduces lip lines and makes lips look fuller and younger looking.

amalian SF 14 natural lips is ideal for cases where mild lip augmentation is needed. It is for age-related volume loss and to remove fine perioral lines.

amalian SF 14 natural lips creates the natural look of young lips without changing the individual's lip geometry.

Needle size: 27G

light lip augmentation





amalian

S & V Technologies GmbH
Neuendorfstraße 20 a
16761 Hennigsdorf
Germany

Tel. +49 (0) 3302-202-7900
Fax +49 (0) 3302-202-7905

info@skin-vision.eu
www.skin-vision.eu